Covid-19 Präventions- konzept 2.0

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
In Vitro Diagnostic Use

Japan's first kit
that can detect
COVID-19 antigens in 30 minutes

SARS-Cov-2 antigen assay kit

ESPLINE[®] SARS-CoV-2

An immunochromatographic assay based on Enzyme Immunoassay (EIA)

<p>1</p>  <p>Antigen</p> <p>Japan's first SARS-CoV-2 antigen assay kit that combines immunochromatographic assay with EIA</p>	<p>2</p>  <p>2 DROPS 1 PRESS</p> <p>Simple procedure: 2 sample drop + 1 press</p>	<p>3</p>  <p>30 min</p> <p>Interpretation time 30 minutes (Less than 30 minutes for 'positive')</p>
--	--	--

Limitations and Precautions

1. Negative results do not preclude SARS-CoV-2 and other virus infection.
2. Do not rely only on the result of this product but also consider clinical symptoms according to the government's guidelines for comprehensive diagnosis.
3. For sample collection and handling, devise necessary measures for biohazards.

SARS-Cov-2 antigen assay kit

ESPLINE® SARS-CoV-2


Read inserted document carefully before use.

Limitations and Precautions

1. Negative results do not preclude SARS-CoV-2 and other virus infection.
2. Do not rely only on the result of this product but also consider clinical symptoms according to the government's guidelines for comprehensive diagnosis.
3. For sample collection and handling, devise necessary measures for biohazards.

SAMPLE COLLECTION AND ASSAY PROCEDURE

SAMPLE COLLECTION

Nasopharyngeal swab fluid specimen :

Insert the sterile swab (Type S) deeply into the nasopharyngeal cavity from nostril to reach adenoids and rub the nasopharyngeal turbinate a few to collect epidermis from the mucosa.

※ To prevent infection, wear a mask, goggles or face shield, proper protective clothing and disposable gloves when collecting samples.

ASSAY PROCEDURE

Collect Samples	Insert the Applicator Tip	Drop Sample	Press	Reaction time
<p>Place the swab into the Sample Extraction Solution (Squeeze Tube) and pinch the swab from outside of the Squeeze Tube approximately 10 times. Finally squeeze the swab head to extract the sample.</p>	<p>Insert the Applicator Tip securely onto the top of the tube of Sample Extraction Solution. <u>Let sample stand for about 5 minutes.</u></p>	<p>Hold Squeeze Tube vertically (approximately 10 mm above the Reaction Cassette) and apply two drops of Squeeze Tube of the sample onto the sample window of the Reaction Cassette.</p>	<p>Press the convex button of the Reaction Cassette down immediately after applying the sample to start the reaction. Leave the cassette horizontally for 30 minutes. Once the reaction is started, keep the interpretation window away from direct wind.</p>	<p>Interpret the result at 30 minutes. If the reference line and the test line appear before 30 minutes, the sample must be considered "positive".</p>

RESULTS

Positive: Both reference line and test line are observed.

Negative: Only reference line is observed.

Retest: ① Only test line is observed. ② Pre-reaction status (red line).

- Interpret the result at 30 minutes. If the reference line and the test line appear before 30 minutes, the sample must be considered "positive". A sample that is "negative" after the 30-minute reaction time and then turns "positive" after 30 minutes, must be considered "negative".
- Retest (1): In case of no blue reference line at 30 minutes or still red reference line at 30 minutes, review the test procedure and repeat the test using a new Reaction Cassette.
- Negative results do not preclude SARS-CoV-2 and other virus infection.

DATA

1 Cross reactivity

No cross reactivity with Influenza viruses as follows;

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No cross reactivity with In-house recombinant human coronavirus antigens as follows;

- MERS-CoV
- HCoV-229E
- HCoV-OC43
- HCoV-NL63
- HCoV-HKU1

2 Comparison with RT-PCR

The test result of this product versus RT-PCR using administrative testing samples is as follows.

		RT-PCR			
		Positive	Negative	Sum	
This product	Positive	16	0	16	Positive concordance rate 66.7% (16/24)
	Negative	8	100	108	Negative concordance rate 100.0% (100/100)
	Sum	24	100	124	Overall concordance rate 93.5% (116/124)

When comparison was based on the number of viral RNA copy numbers, the concordance rates between the two assays was 100% (12/12 cases) for samples with over 1,600 copies/test, 92% (14/15 cases) for samples with over 400 copies or more/test and 83% (15/18 cases) for samples with 100 copies/test or more. Noteworthy, for these samples, the same sample solution was used for the two assays (The swab has been immersed in the virus transport solution then used on ESPLINE SARS-CoV-2)

PRODUCT INFORMATION

Intended Use

for the detection of SARS-CoV-2 antigen directly from nasopharyngeal swab fluid (aid the diagnosis of SARS-CoV-2 infection.)

Product information

CODE NO.	Name	Specimen	Specifications	Expiration	Storage condition
260326	ESPLINE® SARS-CoV-2	Nasopharyngeal swab fluid	100 tests	9 months after	1-30°C

Material provided

Kits components	100 tests
Reaction cassette	10 tests / Box x 10
Sample Extraction Solution(Squeeze Tube)	5 tubes / Bag x 20
< Accessories >Applicator tip	10 tips / Bag x 10
< Accessories >Sterile swab Type S	10 / Bag x 10

[Contact]Global Sales Department
Email : fri_info_s@hugp.com

※ The product specifications are subject to change without prior notice.
※ ESPLINE is a registered trademark of Fujirebio Inc. in Japan and in other countries.

ESL001 2020.6.4000.R-3
(Revised in June, 2020)

www.daswerk.org

www.viennaclubcommission.at